

Dosecare⁺

Hygienic • Customisable • Safe Solution

Unique three-step safety design ensuring the exact dosage time after time.

Find our product launching page at:
www.althenamedical.it




Dosecare⁺ assured fail-safe dosing.

Why Dosecare⁺

- **Fully hygienic system**
Assures the dosing system remains clean.
- **Secure dosing system**
Fail-safe fixed volume allows others to dose safely.
- **Safe and secure attachment**
Safe transfer from adapter to patient.
- **Safe tip design**
Doesn't cause harm to patients mouth.

Version_0.1



Fully customisable heads

- Interchangeable sizes
- Customisable forms
- Full selection of colours

Re-fixable safety switch function

- Possible to re-set the volume
- Impossible to make dosing errors
- Usable without having to check and control every use
- Patent pending system

User-friendly

- Easy and intuitive to utilise
- Easy to clean inside
- Doesn't harm patients mouth

Waste prevention tip


- Fluid doesn't remain on barrel exterior
- Waste prevention technology
- The dose is always what you measure

Adapter bottle attachment

- Fluid doesn't remain inside adapter
- Safe and easy 'click' attachment
- Adapter supplied with kit
- Patent pending system




*Only three steps to fix the correct volume



Other benefits of Dosecare⁺

- User-friendly design
- Fully customisable volumes and colours
- No leakage during dosage phase
- Easy to disassemble and re-assemble
- Doesn't waste content
- Always ensures correct measurement
- 'Click' button for secure dosage
- 'Click' feature for adapter connection
- Blocking design for piston detachment prevention
- Overdose prevention technology

Version_0.2



Every volume setting fixable

- Possible to set every volume, not only a multiple of a fixed one
- Overdose prevention technology
- Interchangeable barrels from v_0.1 to v_0.2
- Patent pending system

Device specifications

Volume*	2.5mL	3mL	4mL	5mL	6mL	7mL	8mL	9mL	10mL
Version	v_0.1					v_0.2			
Parts	Materials								
Piston	Polystyrene					Polypropylene			ABS
Barrel	Polyethylene					Polypropylene			
Adapter	Polyethylene					Polypropylene			
Cursor	Polystyrene					Polypropylene			ABS

All colours available* All volumes available to 12.5mL and more* Branded devices available*

Dosecare⁺ has our company philosophy contained; all devices are highly customisable, in volume, colour and shape. Our quality control is stringent. Every device is strictly checked to comply with pharmaceutical standards. Each device represents who we are and what we stand for.