


Misyri

Markets smallest • Single-use syringe

The smallest syringe on the market. Ideal to be hidden in one hand.

Misyri is engineered to be the smallest syringe on the market. Its' compact and robust design ensures the delivery of medicine to patients, the fastest, shortest and least threatening way possible.


Find our product launching page at:
www.althenamedical.it



Misyri the smallest syringe on the market.

Why Misyri

- **Purposely small**
Designed to be hidden in one hand.
- **Easy vaccinations**
Allowing children to be vaccinated without noticing.
- **Compatible and portable**
The small design reduces packaging space, reducing product and plastic waste.



Fully customisable head

- Interchangeable sizes
- Customisable forms
- Full selection of colours

Removable finger gripper

- Saves transport space
- Customisable form

Fully customisable barrel

- Available in different materials
- Available in different sizes

Luer lock connector

Misyri has our company philosophy contained: the device is highly customisable, in volume, colours and shape. Moreover, quality control is stringent. Every device is strictly checked to comply with pharmaceutical standards. Because each device we produce represents who we are and what we stand for.

Device specifications

Volume	0.2mL	0.3mL	0.4mL	0.5mL	0.6mL
Parts	Materials				
Plunger rod	Polystyrene		Polypropylene		
Barrel	Makrolon®		COP		
Finger gripper	Polyethylene		Polypropylene		
Plunger	Datwyler				

All colours available* All volumes available on request* Branded devices available*