

Sycare

Durable resistant design • innovative closure

Robust quality innovative design. Reliable innovative closure system.

Sycare is engineered to be the strongest and most accurate syringe to date. Its' compact and robust design ensures product protection against air contact.

Find our product launching page at:
www.althenamedical.it



Sycare the most robust syringe on the market.

Why Sycare

- **Highly resistant**
Crack resistant, reinforced screw thread design.
- **Airtight vacuum screw cap**
Engineered air-tight system with friction reduction.
- **Finger grip**
No rotation or movement unless intended.
- **Durable plunger rod**
Doesn't flex allowing high precision and gradual pressure.

Plunger rod

- Durable non-flexible design
- Fully customisable heads

Finger grip


- Supports syringe from movement
- Every form available

Tip design

- Easy-grip design
- Reinforced for crack prevention
- Strong secure screw

Vacuum screw cap

- Easy screw on and off
- No friction between rubber and syringe tip
- Safe and secure when shaken or dropped
- Patent pending system



Why COP

- Transparent like glass
- Weighs less than glass
- Wide range of design options
- Break and leak-resistant
- Cost-effective option
- No alkali ions released
- Consistent PH values

Sycare has our company philosophy contained: the device is highly customisable, in volume, colours and shape. Moreover, quality control is stringent. Every device is strictly checked to comply with pharmaceutical standards. Because each device we produce represents who we are and what we stand for.

Device specifications

Volume*	0.5mL	0.4mL	0.5mL	0.6mL	0.7mL	0.8mL	0.9mL	1.0mL
Parts	Materials							
Plunger rod*	Makrolon®							
Barrel	COP							
Finger grip	Makrolon®							
Screw cap	Makrolon®							
Plunger	Datwyler							

All colours available* All volumes available until 1.0mL and more* Branded plunger heads available*